

# Foreign Body Retrievers

Single Use



 **Medi-Globe**  
The Spirit of Care.™

# Foreign Body Retrievers

Single Use

## Features:

- Wide Opening of Graspers
- Single Use

## Benefits:

- Large Variety of Applications
- Safety and Convenience


## Order Information:

Catalog No.	Description	Outer Diameter	Length
<b>GFB-01-25-180</b>	Foreign Body Retrievers, 4-prong grasper with hooks	2.5 mm	180 cm
<b>GFB-01-25-230</b>	Foreign Body Retrievers, 4-prong grasper with hooks	2.5 mm	230 cm
<b>GFB-02-25-180</b>	Foreign Body Retrievers, 4-prong grasper with eyelets	2.5 mm	180 cm
<b>GFB-02-25-230</b>	Foreign Body Retrievers, 4-prong grasper with eyelets	2.5 mm	230 cm
<b>GFB-03-25-180</b>	Foreign Body Retrievers, forked jaws (2:2) grasper	2.5 mm	180 cm
<b>GFB-03-25-230</b>	Foreign Body Retrievers, forked jaws (2:2) grasper	2.5 mm	230 cm
<b>GFB-04-25-180</b>	Foreign Body Retrievers, forked jaws (1:2) grasper	2.5 mm	180 cm
<b>GFB-04-25-230</b>	Foreign Body Retrievers, forked jaws (1:2) grasper	2.5 mm	230 cm

## Shapes of Foreign Body Retrievers and Polyp Retrievers



Foreign Body Retriever  
4-prong grasper with hooks



Polyp Retriever  
4-prong grasper with eyelets



Foreign Body Retriever  
forked jaws (2:2) grasper



Foreign Body Retriever  
forked jaws (1:2) grasper

Manufactured by:



Medi-Globe GmbH  
Medi-Globe-Strasse 1-5  
83101 Achenmühle, Germany  
Phone: ++49 8032 973 379  
Fax: ++49 8032 973 339  
e-mail: sales@medi-globe.de  
www.medi-globe.de

Distributed by: