


# Biopsy Forceps

Single Use



# Biopsy Forceps

Single Use

## Features:

- Short Rigid Forceps Tips
- Large Jaw Opening
- Single Use

## Benefits:

- Easy Passage through Angled Working Channels
- Facilitates Tangential Biopsies
- Safety and Convenience

## Order Information:

Catalog No.	Description	Outer Diameter	Length	Minimum Working Channel Diameter
<b>GBF-01-22-180</b>	Biopsy Forceps, oval cups	2.4 mm	180 cm	2.6 mm
<b>GBF-01-22-230</b>	Biopsy Forceps, oval cups	2.4 mm	230 cm	2.6 mm

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH  
Medi-Globe-Strasse 1-5  
83101 Achenmühle, Germany  
Phone: ++49 8032 973 379  
Fax: ++49 8032 973 339  
e-mail: sales@medi-globe.de  
www.medi-globe.de

Distributed by: